

➤ MEDICAL IMAGING SYSTEM

Digital X-RAY Imaging System CJ-DR0



Product Features

- Support linkage with hospital management software
- Support consultation room and DR room interoperability
- Support for consultation room film reading
- Complete software tools
- Clear images
- Easy operation with dual screen display (optional)



➤ Product Main Parts



**Microfocus bulb tube (optional Canon bulb tube)**  
Difficult images at a glance, suitable for cat specialist hospitals, reducing artifacts in animals due to shortness of breath, heartbeat, nervous tremors, etc., avoiding misdiagnosis and missed diagnosis.



**17\*17 fixed flat detector**  
Advanced cesium iodide evaporation technology, fast imaging speed, multiple trigger modes, 16-bit conversion, 3.6Lp/mmg high spatial resolution, better image quality and lower exposure dose



**X-Ray Generator**  
32KW high power, nominal d current up to 400ma, large pets, small pets perfect for  
Dynamic and static high precision voltage and current replenishment

➤ PATIENT MONITORING & LIFE SUPPORT


Infusion Pump CJ-HF710B



➤ Product Features

- Adjustable Levels Setting:10 levels settingfor occlusion, air bubble and alarm volume.
- Strong Battery Backup:3000mAh battery and backup time more than 10hrs.
- User Friendly: Big Screen and numeric keypad.
- Dual Occlusion Sensor: Up and down occlusion sensor with DPS system.
- Infusion Mode:Volume mode,drip mode,time mode.

➤ Application Scenarios



| Item                        | Specification                                                                                                      | Item                        | Specification                                                                                                                                                                                                                                                                                                       |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Infusion Mode               | Volume Mode, Drip Mode, Time Mode, Drug Mode                                                                       | Air Bubble Detector         | Ultrasonic detection, ≥50ul , Single bubble volume: 50~300 ul,10 levels adjustable                                                                                                                                                                                                                                  |
| Rate Unit                   | ml/h, volume/time, d/min                                                                                           | Volume                      | 10 levels adjustable                                                                                                                                                                                                                                                                                                |
| IV Sets                     | Standard single use IV sets (10,15,20,60d/mL, with diameter3.4~4.5mm)                                              | Motor Anti-reverse Function | Anti-reverse function to prevent upstream                                                                                                                                                                                                                                                                           |
| Flow Rate                   | 1.0ml/h~1200ml/h; 1.0d/min~400d/min                                                                                | IV Set Management           | Open system, accept different IV brands after calibration                                                                                                                                                                                                                                                           |
| Increment                   | 1.0-99.9ml/h in 0.1ml/h step<br>100-1200ml/h in 1ml/h step<br>1-400d/min in 1d/min step                            | Self-Testing                | Self-Testing before startup and during infusion                                                                                                                                                                                                                                                                     |
| Deviation                   | Flow Rate: ≤5% ( after calibration)                                                                                | Audiovisual Alarm Message   | Air bubble, upstream occlusion, downstream occlusion, door open, infusion finish, standby time expired, battery low, battery empty, battery disengaged, AC Power disengaged, motor error, pump error; pressure sensor error, parameter error; drip sensor error, drip sensor disengaged, bottle empty, system error |
| VTBI                        | 0.1 ~9999.9 ml                                                                                                     | Screen                      | 3.2inch brightness LCD                                                                                                                                                                                                                                                                                              |
| Total Volume                | 0.0~9999.9 ml                                                                                                      | Power Supply                | AC Power : 100V~240V, 50Hz/60Hz, ≤25VA                                                                                                                                                                                                                                                                              |
| Time                        | 1min~99h00min                                                                                                      | Battery                     | Rechargeable Lithium battery, DC11.1V, 3000mAh operating time: ≥10h @25ml/h, displaying remaining battery                                                                                                                                                                                                           |
| Purge Rate                  | 600ml/h, display purge volume                                                                                      | Operation Condition         | Temperature: +5°C~+40°C, Humidity : 20% ~90%, Air pressure : 70~106KPa                                                                                                                                                                                                                                              |
| KVO                         | 1.0ml/h~5.0ml/h, adjustable                                                                                        | Transport & Storage         | Temperature :-20°C~+55°C,Humidity: 10%~95%, Air pressure:50~106KPa                                                                                                                                                                                                                                                  |
| Keypad                      | Numeric keypad                                                                                                     | Classification              | Type CF , Class I ,IPX3                                                                                                                                                                                                                                                                                             |
| BOLUS                       | Bolus function available, flow rate adjustable                                                                     | Dimension                   | 145mm X 150mm X 200mm; Weight: 1.7kg                                                                                                                                                                                                                                                                                |
| Night Mode                  | Screen to be dim automatically in 5mins after stand by time expired                                                |                             |                                                                                                                                                                                                                                                                                                                     |
| Occlusion Pressure Detector | Upstream and downstream pressure detection, 40KPa~130KPa, 10 levels adjustable with dynamic pressure value display |                             |                                                                                                                                                                                                                                                                                                                     |

Veterinary

IVD Diagnostic Solution

Medical Imaging System

Patient Monitoring & Life Support



Healgen Scientific Limited Liability Company

Tel: +1 713-733-8088 Email: sales.healgen01@healgen.com Website: www.healgen.com

Address: 3818 Fuqua Street, Houston, TX 77047, USA.

 CJ-G200 Vet fluorescence immunoassay analyzer



 Instant printer

 On-Site detection

 Portable

 Touch screen operation

 ID chip

 Sample testing

Fast ≤8s

 Application Scenarios



**Veterinary Institution**

- Physical Examination Center
- Veterinary School
- Pet Shop



**Animal Hospital**

- Outpatient/emergency labs
- Clinical Laboratory
- Laboratory
- Ambulance

 Technical Parameters

| Item         | Specification                                     | Item               | Specification                                            |
|--------------|---------------------------------------------------|--------------------|----------------------------------------------------------|
| Power Supply | 12V, 5A                                           | Network            | Support RS232 serial port, wireless Wi-Fi, Bluetooth, 4G |
| Size         | 232mm×269mm×147mm                                 | Detection speed    | 8 seconds                                                |
| Weight       | 2050g                                             | Storage capacity   | 10000 groups datas                                       |
| Display      | 7-inch LCD capacitive screen, 1024×600 resolution | Detection method   | Support one card with three items                        |
| Printer      | Thermal printer, support thermal induction        | Language interface | Chinese and English, other languages can be expanded     |

 Fluorescence Rapid Test Kit

|                                                                                       | Category              | Product Name                                   | PN.     | Sample                          |
|---------------------------------------------------------------------------------------|-----------------------|------------------------------------------------|---------|---------------------------------|
|    | Inflammatory Diseases | cCRP Rapid Test Kit ( Fluorescence )           | CJ-F001 | Serum/Plasma/Whole Blood        |
|                                                                                       | Intestinal Diseases   | CPV Ag Rapid Test Kit ( Fluorescence )         | CJ-F003 | Feces                           |
|                                                                                       |                       | CCV Ag Rapid Test Kit ( Fluorescence )         | CJ-F004 | Feces                           |
|                                                                                       |                       | CCV/CPV Ag Rapid Test Kit ( Fluorescence )     | CJ-F005 | Feces                           |
|                                                                                       |                       | CPL Rapid Test Kit ( Fluorescence )            | CJ-F002 | Serum/ Plasma                   |
|                                                                                       | Hormones              | cProg Rapid Test Kit ( Fluorescence )          | CJ-F006 | Serum/Plasma                    |
|                                                                                       |                       | cT4 Rapid Test Kit ( Fluorescence )            | CJ-F021 | Serum/Plasma                    |
|                                                                                       |                       | cTSH Rapid Test Kit ( Fluorescence )           | CJ-F025 | Serum/Plasma                    |
|                                                                                       |                       | cCor Rapid Test Kit ( Fluorescence)*           | CJ-F026 | Serum/Plasma                    |
|                                                                                       | Respiratory Diseases  | CDV Ag Rapid Test Kit ( Fluorescence )         | CJ-F008 | Secretions (Eyelids/Mouth/Nose) |
|                                                                                       |                       | CAV Ag Rapid Test Kit ( Fluorescence )         | CJ-F007 | Secretions(Mouth/Nose)          |
|                                                                                       | Cardiac Diseases      | cD-Dimer Rapid Test Kit ( Fluorescence)*       | CJ-F027 | Plasma                          |
|                                                                                       | Kidney Diseases       | cSDMA Rapid Test Kit ( Fluorescence)           | CJ-F023 | Serum/ Plasma                   |
|  | Others                | CRV Ab Rapid Test Kit ( Fluorescence )         | CJ-F009 | Serum/ Plasma                   |
|                                                                                       |                       | CRV Ag Rapid Test Kit ( Fluorescence )         | CJ-F010 | Secretions(Mouth/Nose)          |
|                                                                                       |                       | CDV/CPV/ICH Ab Rapid Test Kit ( Fluorescence ) | CJ-F011 | Serum/ Plasma                   |
|                                                                                       | Inflammatory Diseases | fSAA Rapid Test Kit ( Fluorescence )           | CJ-F012 | Serum/Plasma /Whole Blood       |
|                                                                                       | Intestinal Diseases   | FPV Ag Rapid Test Kit ( Fluorescence )         | CJ-F014 | Feces                           |
|                                                                                       |                       | FPL Rapid Test Kit ( Fluorescence )            | CJ-F013 | Serum/ Plasma                   |
|                                                                                       |                       | FCoV Ag Rapid Test Kit ( Fluorescence )        | CJ-F015 | Serum/ Plasma                   |
|                                                                                       | Hormones              | fT4 Rapid Test Kit ( Fluorescence )            | CJ-F022 | Serum/ Plasma                   |
|                                                                                       | Respiratory Diseases  | FHV Ag Rapid Test Kit ( Fluorescence )         | CJ-F016 | Secretions (Mouth/Nose)         |
|                                                                                       |                       | FCV Ag Rapid Test Kit ( Fluorescence )         | CJ-F017 | Secretions (Mouth/Nose)         |
|                                                                                       | Cardiac Diseases      | fNT-proBNP Rapid Test Kit ( Fluorescence )     | CJ-F018 | Serum/ Plasma                   |
|                                                                                       | Kidney Diseases       | FSDMA Rapid Test Kit ( Fluorescence )          | CJ-F024 | Serum/ Plasma                   |

Products marked with "\*" will be coming soon

➤ IVD SOLUTION

 Fluorescence Rapid Test Kit

|                                                                                     | Category | Product Name                                   | PN.     | Sample        |
|-------------------------------------------------------------------------------------|----------|------------------------------------------------|---------|---------------|
|  | Others   | FRV Ab Rapid Test Kit ( Fluorescence )         | CJ-F019 | Serum/ Plasma |
|                                                                                     |          | FPV/FHV/FCV Ab Rapid Test Kit ( Fluorescence ) | CJ-F020 | Serum/ Plasma |

 Colloidal Gold Rapid Test Kit

|                                                                                       | Category             | Product Name                                               | PN.     | Sample                          |
|---------------------------------------------------------------------------------------|----------------------|------------------------------------------------------------|---------|---------------------------------|
|    | Intestinal Diseases  | CPV Ag Rapid Test Kit ( Colloidal gold )                   | CJ-C001 | Feces                           |
|                                                                                       |                      | CCV Ag Rapid Test Kit ( Colloidal gold )                   | CJ-C002 | Feces                           |
|                                                                                       |                      | CCV/CPV Ag Single Window Rapid Test Kit ( Colloidal gold ) | CJ-C003 | Feces                           |
|                                                                                       |                      | CCV/CPV Ag Double Window Rapid Test Kit ( Colloidal gold ) | CJ-C004 | Feces                           |
|                                                                                       | Respiratory Diseases | CDV Ag Rapid Test Kit ( Colloidal gold )                   | CJ-C005 | Secretions (Eyelids/Mouth/Nose) |
|                                                                                       |                      | CIV Ag Rapid Test Kit ( Colloidal gold )                   | CJ-C006 | Secretions (Mouth/Nose)         |
|                                                                                       |                      | CAV Ag Rapid Test Kit ( Colloidal gold)                    | CJ-C007 | Secretions (Eyelids/Mouth/Nose) |
|                                                                                       |                      | CDV/CAV Ag Rapid Test Kit ( Colloidal gold)                | CJ-C008 | Secretions (Eyelids/Mouth/Nose) |
|                                                                                       |                      | CIV/CDV/CAV Ag Rapid Test Kit ( Colloidal gold )           | CJ-C009 | Secretions (Eyelids/Mouth/Nose) |
|  | Intestinal Diseases  | FPV Ag Rapid Test Kit ( Colloidal gold )                   | CJ-C010 | Feces                           |
|                                                                                       |                      | FCoV Ag Rapid Test Kit                                     | CJ-C011 | Feces                           |
|                                                                                       |                      | FPV/FCoV Ag Rapid Test Kit                                 | CJ-C012 | Feces                           |
|                                                                                       | Respiratory Diseases | FHV Ag Rapid Test Kit ( Colloidal gold )                   | CJ-C013 | Secretions (Mouth/Nose)         |
|                                                                                       |                      | FCV Ag Rapid Test Kit ( Colloidal gold )                   | CJ-C014 | Secretions (Mouth/Nose)         |
|  | Parasitic Diseases   | TOXO Ab Rapid Test Kit                                     | CJ-C015 | Serum/ Plasma                   |